

Stericlin® Instrument bags



protect and fix your instruments

Forceps, knife handles, probes, fine accessories etc. are often added slackly to the instrument tray inside a kidney basin. Thus high-quality instruments may easily get damaged during transport. Moreover, components which have fallen out of the kidney basin are tediously to find or to remove out of the tray.

Stericlin® Instrument bags protect your products absolutely reliably from damages by the combination of sterilisation-resistant non-woven and transparent laminate film. This material-combination has proved itself millions of times with sterile barrier systems.

The division of the bag into separate parts and different bag sizes safeguards a clear arrangement, quick identification and intentional removal of the assorted instruments.

Attention! This product is not a sterile barrier system! Upon assembling the instruments you can fold or roll the bags, add them to the instrument tray and sterilise them in the container (or in flexible or see-through-packaging). Instrument bags are intended for single-use only.



Instrument bags

Article-No.	Size in cm	Pack (pcs.)
3FSZB730120	52 x 21	250